

IMIDASECT ANTS



PROFESSIONAL USE

PROFESSIONAL USE

PROFESSIONAL USE

COMPOSITION

Imidacloprid 0,01% as gel

DESCRIPTION

IMIDASECT ANTS is an insecticide bait formulated as a gel, to control specific feeding ants sugars and come to constitute a pest, as Argentina ant (*Linepithema humile*), and those considered harmful as *Monomorium pharaonis*, *Tapinoma Melanocephalum*, *Lasiud Niger*, *Lasius flavus*, *Camponotus pennsylvanicus*.

PLACE OF APPLICATION AND PESTS

For use indoors and outdoors. Timely implementation of the product in cracks and crevices in the transit zones and / or proliferation of insects, by syringe as drops. IMIDASECT ANTS are designed for the controlled placement of bait in the following situations: Pharaoh Ants (*Monomorium pharaonis*); Argentine Ants (*Linepithema humile*); Black Ants (*Lasius niger*).

Domestic premises (including kitchens)

Food Handling: Processing: Food manufacture, kitchens, slaughterhouses etc. Storage: Food retailers, warehouses, raw material stores, silos (except where grain is stored) etc. Preparation: Restaurants, commercial kitchens, etc.

Public Buildings: Continuously occupied areas including hospital wards, hotels, public baths, municipal buildings, churches, halls, community centers, cinemas, etc.

MODE OF APPLICATION

Direct Application from the commercial pack, application in localized ways in way of small drops

DOSES

Indoor Dose: Apply 1-2 drops of gel per m², 1 drop of gel per meter for black ants.

The droplets should be approximately 1 cm (equivalent to 0.2 g)

Outdoor Dose: Apply 1 drop of gel directly at the entrance of the nests or ant trails transit.

The droplets should be approximately 1 cm (equivalent to 0.2 g)

OBSERVATIONS

The drops will remain actives for several months, although it is advisable to check after 1 week depending the degree of infestation to recharge in case the bait has been consumed. Can be used both inside (Kitchen, Stockrooms, warehouse, so on) and outside (patio, terrace, balcony) . In those places where proliferate or move the ants.

PRESENTATION

Size	Presentation	Showing	Units x case	Cases x Box
35 gr.	Syringe of HDPE transparent without pushed	Individual labeled syringe	6 (Likely could be also 12)	6